

SENSOVISC

*Viscoelastic Gel
for Intra-Articular Injection*

✓ Sterile

✓ Viscoelastic

✓ Non-pyrogenic

✓ Highly purified
sodium hyaluronate


Istem

MEDİKAL TIBBİ CİHAZ VE SANAYİ LTD.ŞTİ.

CE 2292

Medical Devices Class III


SENSOVISC

Sensovisc is a synovial fluid substitute which, thanks to its viscoelastic and lubricant properties, promotes the restoration of rheological conditions of the joints, altered in degenerative or post-traumatic conditions. The sterile sodium hyaluronate gel, improving the characteristics of the synovial fluid, exerts a protective action of the joints and helps the improvement of joint function and the reduction of pain symptoms. **Sensovisc** acts only at the joint where it is injected, without exerting any systemic action.



Sodium hyaluronate

- ✓ High molecular weight 2.1-2.4 million daltons
- ✓ A biodegradable polymer
- ✓ Obtained by bacterial fermentation
- ✓ No animal origin



SENSOVISC



SENSOVISC 75mg

25 MG HA(2.5%) IN 3 ML PHOSPHATE BUFFER SOLUTION

- High Molecular Weight Sodium Hyaluronate
- Bacterial Fermentation
- pH 6.8-7.8
- Osmolality 270-400 mOsm/kg, similar to the synovial fluid
- Viscosity 200-250.000 cP

SENSOVISC 60mg

20 MG HA(2.0%) IN 3 ML PHOSPHATE BUFFER SOLUTION

- High Molecular Weight Sodium Hyaluronate
- Bacterial Fermentation
- pH 6.8-7.8
- Osmolality 270-400 mOsm/kg, similar to the synovial fluid
- Viscosity 100-150.000 cP



SENSOVISC 40mg

20 MG HA(2.0%) IN 2 ML PHOSPHATE BUFFER SOLUTION

- High Molecular Weight Sodium Hyaluronate
- Bacterial Fermentation
- pH 6.8-7.8
- Osmolality 270-400 mOsm/kg, similar to the synovial fluid
- Viscosity 100-150.000 cP

SENSOVISC 32mg

16 MG HA(1.6%) IN 2 ML PHOSPHATE BUFFER SOLUTION

- High Molecular Weight Sodium Hyaluronate
- Bacterial Fermentation
- pH 6.8-7.8
- Osmolality 270-400 mOsm/kg, similar to the synovial fluid
- Viscosity 40-120.000 cP



SENSOVISC 20mg

10 MG HA(1.0%) IN 2 ML PHOSPHATE BUFFER SOLUTION

- High Molecular Weight Sodium Hyaluronate
- Bacterial Fermentation
- pH 6.8-7.8
- Osmolality 270-400 mOsm/kg, similar to the synovial fluid
- Viscosity 20-60.000 cP

SENSOVISC

| REFERENCE NO | PRODUCT NAME | DEFINITION |
|--------------|-----------------|------------------------------------------------------------------|
| SV-2-20 | SENSOVISC 20 MG | Sterile Sodium Hyaluronate Gel- 2 mL Prefilled Syringe (1.0%) |
| SV-2-32 | SENSOVISC 32 MG | Sterile Sodium Hyaluronate Gel- 2 mL Prefilled Syringe (1.6%) |
| SV-2-40 | SENSOVISC 40 MG | Sterile Sodium Hyaluronate Gel- 2 mL Prefilled Syringe (2.0%) |
| SV-3-60 | SENSOVISC 60 MG | Sterile Sodium Hyaluronate Gel- 3 mL Prefilled Syringe (2.0%) |
| SV-3-75 | SENSOVISC 75 MG | Sterile Sodium Hyaluronate Gel- 3 mL Prefilled Syringe (2.5%) |